# Digital microscope operating instructions

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
1. USB接口(The USB interface) 2.LED灯光(LED lights) 3.拍照按钮(Photo button)  
4.调控灯光亮度(Adjust light brightness) 5.U型调节支架(U-shaped adjusting support)  
6.调焦滚轴(Focusing roller) 7.高度、角度调节旋钮(Height and Angle adjusting knob)


# Measuring software

## Windows installation software method

1. Search and click <http://soft.hvscam.com>. Please select the version you need windows/mac/ acmcap. Then download the software, as shown in the figure:

2. Double click and open. The following information will appear. Then, click next step:

3. Reach the following license agreement, click "accept":

# Measuring software

## Windows installation software method

4. Reach the following information select installation site and click install:

5. Reach the information of being installed, wait a minute, wait for automated installer:


6. Reach the following information: select and operate HView 1.0 (R), then click "finish", PC automatically run a software.

# Using the software

## The microscope comes into use

1. Connect your microscope to PC via the USB cable.  open the software you can see the flowing interface after stepping into the software:

2. Click "equipment" in the menu bar. Then, click "open", the following information will appear:


3. As the picture shown above, select the microscope name "GENERAL-UVC". Then select "color/compression" as the local computer "MJPG". According to the self-adaption of microscope MJPG, YUY2 or select manually, select the image, after selecting the image resolution, click "open" into the preview screen, as shown in the figure below:

## Measuring software

### Mac installation software method

1. Search and click <http://soft.hvscam.com>. Please select the version you need windows / mac / acmcap. Then download the software, as shown in the figure:

2. Double click and open HiViewSetu.dmg. The following information will appear. Then, click next step:

## Using the software

### The microscopes comes into use

1. Connect your microscope to a mac computer or laptop via the USB cable. Then, double click the desktop  open the software. If there are safety tips:

# The solution is as follows

## 1. System Preference

## 2. Open the Security & Privacy Select App Store and admmissive developer

Select "Still to open"


3. The opening interface of software is as shown in the figure.

Click  in the menu, open the image Reach HiView If you want to access the camera Select Yes

# Menu bar function introduction


## Measurement:

Click to perform operations including cross curve, straight line, parallel line, circle, concentric circle, rectangle, Angle, ellipse, arc, polygon, arbitrary drawing, words, stop measurement, delete recent, delete all, calibration, and options under the preview screen and photo screen.


 Calibration -> Modulate the microscope image into a clear picture, and place the correction bar at the extremitas anterior of the microscope for correction, as shown in the figura below.  For example, if the specification of calibration bar is 0.1mm per small grid, the software calibration bar should be in the 22nd small grid of the calibration bar and input 2.2.

The calculation formula is  $0.1\text{mm} \times 22\text{mm} = 2.2\text{mm}$

(Note: the object distance shall be recalibrated for each change, and the measurement of each item can be carried out after the recalibration.)

Options-> The line color, line width, word color and word size can be set, as shown in the figure:

## Assistance:

Language → Click to switch among simplified Chinese, English, Japanese, Russian, French, German, Italian, Spanish, Portuguese, Korean and other languages  
About (A) → Click to view the software information and version number.

## The Function Declaration of Shortcut Key

 Open device options

 video recording (video)

 1280x720 modify the resolution ratio

 show cross line of reference

 draw parallel lines

 draw concentric circle

 draw angle

 draw arc

 arbitrary drawing

 stop measurement

 return to previous operation

 Set the line width, the size and color of word

 Grab images (photograph) (Support space bar shot out photo)

This software does not support hardware photos. If you need hardware photos, please download Arcaap v3.0.0

 100% modify the display scale

 show cross line of reference

 draw circle

 draw rectangle

 draw ellipse

 draw polygon

 mark words

 microscope calibration

 delete all measurement

### Photo Files

Image	Video
01/05/2019	
09:31 001	
09:38 001	
09:38 002	
09:38 003	
11:20 001	

### Audio Record Files

Image	Video
01/05/2019	
11:25 002	
11:26 001	

If there is any question, please contact our company. The company reserves the right of final interpretation.

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