
Levenhuk DTX 800 LCD

Digital Microscope

EN	User Manual
BG	Ръководство за потребителя
CZ	Návod k použití
DE	Bedienungsanleitung
ES	Guía del usuario
HU	Használati útmutató
IT	Guida all'utilizzo
PL	Instrukcja obsługi
PT	Manual do usuário
RU	Инструкция по эксплуатации
TR	Kullanım kılavuzu

levenhuk 
Zoom&Joy


Levenhuk DTX 800 LCD

Digital Microscope

EN

Microscope parts

- | | |
|--|---|
| 1 MicroSD card slot | 10 Digital zoom |
| 2 HDMI out | 11 Brightness adjustment wheel (lower illumination) |
| 3 Focusing wheel | 12 5V/USB Type-C port |
| 4 Polarizer wheel | 13 Power on/off |
| 5 Charging indicator | 14 Settings |
| 6 Lower illumination | 15 Up |
| 7 Height adjustment knobs | 16 Down |
| 8 Capture | 17 OK |
| 9 Brightness adjustment wheel (upper illumination) | 18 Photo/Video/Playback |

The kit includes: microscope, USB Type-C cable, HDMI cable, power adapter, software CD, calibration slide, user manual, and lifetime warranty.

Never, under any circumstances, look directly at the Sun, another bright source of light or at a laser through this device, as this may cause **PERMANENT RETINAL DAMAGE** and may lead to **BLINDNESS**.

Microscope assembly

The battery has to be fully charged before turning the microscope on for the first time. Plug the USB Type-C cable into the power cable socket on the microscope body (12) and connect it via the power adapter to a power source. The charging indicator (5) will turn red. Once the battery is fully charged, the indicator will go out. You can now disconnect the microscope from the power source and turn it on. The fully-charged battery allows you to use the microscope for two hours. The illumination will grow dim and the image quality will decrease to inform you that the battery is low. Charge the battery as described above.

A microSD card (purchased separately) is required to save captured images and recorded videos. Turn the microscope off and insert the microSD card into the corresponding slot on the microscope body (1) until it locks in place with a click. Do not apply excessive force when inserting the microSD card. Try rotating the card if it does not lock in place with ease. Format the card before using it. To turn the microscope on, press the power button (13). To turn the microscope off, press the power button (13) again (refer to «Power saving» for more information).

Display icons

Capture an image

Playback mode

12M

Current image resolution

Battery status indicator

Start and stop video recording

1X

Magnification, x

MicroSD status icon

Capturing an image

Turn the microscope on and place a specimen on the stage. The height of the objective lens can be regulated with height adjustment knobs (7). Adjust the illumination brightness, so that the specimen is evenly lit. You can use upper or lower illumination. Change magnification and focus your view (3), if necessary. Press the Capture button (8).

Recording a video

Press the Photo/Video/Playback button (18) to switch to video mode. The icon in the upper left corner of the screen will change. Press the Capture button (8) to start recording and press it again to stop.

Playback

Press the Photo/Video/Playback button (18) to switch to playback mode. The icon in the upper left corner of the screen will change. Select the desired photo or video clip with Up (15) and Down (16) buttons and press the OK button (17) to play it.

Polarizing filter

You can use a polarizing filter to reduce light reflection while observing objects with a shiny surface like metal, crystals, etc. Rotate the polarizing wheel (4) by 90° to activate or deactivate the polarizing effect.

Mass storage on PC

Insert a microSD card in the slot. Connect the microscope to your PC with the USB Type-C cable to download images from the card to your computer.

Preview on PC

Only works with the installed software (included). Remove the microSD card from the slot. Connect the microscope to your PC with the USB Type-C cable.

Connecting the microscope to a TV set

Plug the HDMI cable into the HDMI out (2) of the microscope and connect it to your TV set (refer to the user guide of your TV set for additional information). Turn on the TV and select the correct channel from the corresponding menu of the TV set (refer to the user guide of your TV set for additional information).

Display menu

Press the Settings button (14) to open the main display menu. In this menu you can set image and video resolution, the interface language, date and time, turn on the power saving mode, format the SD card and reset settings. Also, this menu provides access to the Loop recording mode and Time Stamp option.

Image resolution

Press the Settings button (14) to select **Photo Resolution** from the main menu, then choose the required resolution by pressing the Up (15) or Down (16) buttons, and confirm by pressing the OK button (17).

This model supports the following image resolutions: 16M, 12M, 10M, 8M, 5M, 3M.

Capture mode

Press the Settings button (14) and choose **Time Lapse Photo** to setup the sequence capture. For example, you can set up the microscope to take a photo every 5 seconds.

Video resolution

Press the Settings button (14) to check the **Video Resolution**.

This model supports the video resolution 1440P.

Loop recording

Press the Settings button (14) and choose **Record Cycle**. You can select the duration of the video clip or turn the Loop recording mode off. While the Loop recording mode is on, each video clip of selected duration will be recorded over the previous one.

Time stamp

Press the Settings button (14) and choose **Time Stamp** to select a target option.

Cursor mark switch

Press the Settings button (14) and choose **Reticle Switch** to select a target option.

Date and time

Press the Settings button (14) and choose **Date/Time**, then set the current date and time by pressing the Up (15) or Down (16) buttons, and press (18) to confirm settings.

Language

Press the Settings button (14) and choose **Language**, then select one of the existing settings. Press the OK button (17) to confirm selection.

Power saving

Press the Settings button (14) and choose **Auto Power Off** to select the idle time before the microscope automatically shuts down. Press the OK button (17) to confirm selection.

MicroSD card formatting

Press the Settings button (14) and choose **Format** to format the installed memory card. All the files on the card will be deleted after card formatting. As a rule, formatting the memory card is not necessary, unless your card is not able to save the captured image.

Settings reset

Press the Settings button (14) and choose **Default Setting** to revert back to the factory settings.

PortableCapture Plus software

Insert the installation CD into your CD drive. Choose the folder according to the operating system of your device. Double click the installation file to install the application and follow the instructions on the screen. The **User Manual** folder on the installation CD contains a copy of the user manual. Connect the microscope to your PC with a USB Type-C cable and run the PortableCapture Plus application. A pop-up message will warn you if the microscope is not connected to the PC.

Application menu

Capture an image

Start and stop video recording

Sequence capture (you can set the time to begin capturing images or videos, interval between them, video duration and the number of images or videos to be captured)

Full screen mode (to exit full screen mode, press Esc on your keyboard or double-click anywhere on the screen)

File

Photos Folder: set a directory for captured images.

Videos Folder: set a directory for recorded videos.

Options

Resolution: set image resolution.

Date/Time: hide or display date and time of capture during preview.

Language: change the language of the user interface.

XY Guide: hide or display a grid on images.

Full screen mode: go full screen.

Capture

Photo: capture an image.

Video: record a video.

Capturing an image

Click the corresponding icon on the toolbar or select **Photo** from the **Capture** menu.

Recording a video

Click the corresponding icon on the toolbar or select **Video** from the **Capture** menu.

Sequence capture

Click the corresponding icon on the toolbar. A dialog window will appear where you can set a number of parameters for sequence capture. Click the **Photo** radio button to set the time to begin capturing images, interval between them and the total number of images. Click the **Video** radio button to set the time to begin recording videos, their duration, interval between them and the total number of clips.

Save as

Copy to clipboard

Copy file

Previous image

Next image

Undo the last action

Repeat the last action

Draw

Text

Measurement

Calibration

Calibration

Before calibrating the microscope, turn the grid display on. To do so, select **On** from **Options > XY Guide**. By default, the grid is visible. Place the calibration scale on the stage and focus your view. Make sure that the vertical axis of the grid is parallel to scale divisions, and capture an image. Double-click the image thumbnail to open it in a separate window.

Click the corresponding icon on the toolbar to begin calibration. Using your mouse, choose two points on the image (you have to know the actual distance between these two points). Once the second point is chosen, enter the known distance in the **Actual dimension** field in the dialog window. The application will automatically calculate the magnification of the image. Click OK to close the dialog window. The calculated magnification will be displayed in the **Magnification** field.

Chosen points should form a horizontal line.

You can check calibration results by measuring the same distance with the **Straight line** tool. Click the measurements icon, choose the **Straight line** tool and draw a similar line on the image. If the measured distance is equal to the actual distance, calibration has been successful.

Recalibration

The calibration process has to be repeated if the magnification or focus has been changed during observations.

Click the corresponding icon on the toolbar again and select **Reset picture magnification** from **Reset magnification**. Repeat the calibration process, as described above. Change the magnification and capture a few images. Double-click one of the thumbnails to open the resulting image in a preview window. Click the calibration icon on the toolbar and select **Set picture magnification** from **Set magnification**. Enter the current magnification in the dialog window and click OK.

If you want to use the same magnification when capturing images in the future, click the calibration icon on the toolbar and select **Set capture magnification** from **Set magnification**.

Measurements

It is recommended to calibrate the system before proceeding with measurements.

Straight Line. Press and hold the left mouse button to start drawing a line. Release the left mouse button to complete the line.

Continuous Line. Press and hold the left mouse button and draw a line. Release the left mouse button to complete the segment. You can continue adding more segments to the free-form line.

Radius Circle. Draw a straight line on the image. The application will automatically calculate the radius, as well as circumference and area of the corresponding circle.

Diameter Circle. Draw a straight line on the image. The application will automatically calculate the diameter, as well as circumference and area of the corresponding circle.

Three Points Angle. Press and hold the left mouse button to draw a line on the image. Release the left mouse button and draw another line. The application will automatically calculate value of the resulting angle.

Draw and text

PortableCapture Plus allows you to add drawings and textboxes to your images. Open an image in a preview window and click the **Pencil** icon on the toolbar. Select one of the tools from the drop-down menu and draw something on the image. By clicking the TT icon on the toolbar, you can edit font type and color to be used in textboxes.

Specifications

LCD display	7", 1024x600px
Optics material	optical glass
Megapixels	4
Magnification	20–300x, up to 1500x with digital zoom
Photo	*.jpg
Video	*.mp4
Focusing	manual, 5–70mm
Speed	30fps
Illumination	LEDs, with brightness adjustment
Body material	plastic
Power source	5V DC-in via USB Type-C cable or built-in Li-Ion battery 3.7V, 9.25Wh, 2500mAh; battery duration: 2 hours; charging time: 4.5 hours
Software	image and video processing software with measurement function
Software language	English, German, French, Spanish, Russian, Japanese, Chinese, etc.
Operating system	Windows 10/11, Mac 10.15 and above
Ability to connect additional equipment	support of microSD cards with capacity up to 128GB (not included in the kit); connect to PC via USB Type-C cable (included); connect to TV via HDMI cable (included)
Operating temperature range	0... +40°C / +32... 104°F

The manufacturer reserves the right to make changes to the product range and specifications without prior notice.

Care and maintenance

White LEDs are very bright; do not look at them as this may cause permanent damage to your vision. Take necessary precautions when using the device with children or others who have not read or who do not fully understand these instructions. Do not try to disassemble the device on your own for any reason. For repairs and cleaning of any kind, please contact your local specialized service center. After unpacking your microscope and before using it for the first time check for integrity and durability of every component and connection. Do not touch the optical surfaces with your fingers. To clean the device exterior, use only special cleaning wipes and special optics cleaning tools from Levenhuk. Do not use any corrosive or acetone-based fluids to clean the optics. Abrasive particles, such as sand, should not be wiped off lenses, but instead blown off or brushed away with a soft brush. Do not apply excessive pressure when adjusting focus. Do not overtighten the locking screws. Protect the device from sudden impact and excessive mechanical force. Store the device in a dry, cool place away from hazardous acids and other chemicals, away from heaters, open fire and other sources of high temperatures. Do not use the device for lengthy periods of time, or leave it unattended in direct sunlight. Keep the device away from water and high humidity. Do not submerge. Be careful during your observations, always replace the dust cover after you are finished with observations to protect the device from dust and stains. To avoid data loss or equipment damage, always close the application first, and then disconnect the microscope from your PC. Wait until the illumination goes off completely, and disconnect the microscope. Children should use the device under adult supervision only. **Seek medical advice immediately if a small part or a battery is swallowed.**

Battery safety instructions

The device is equipped with a rechargeable lithium-ion battery. This avoids frequent battery replacement. Always switch the device off when not in use. If the battery charge is low, please recharge the device in time. Do not overheat the battery. Do not discharge the battery completely. Keep batteries out of the reach of children, to avoid risk of ingestion, suffocation, or poisoning. Utilize used batteries as prescribed by your country's laws.

Levenhuk International Lifetime Warranty

All Levenhuk telescopes, microscopes, binoculars, and other optical products, except for their accessories, carry a lifetime warranty against defects in materials and workmanship. A **lifetime warranty** is a guarantee on the lifetime of the product on the market. All Levenhuk accessories are warranted to be free of defects in materials and workmanship for **six months** from the purchase date. The warranty entitles you to the free repair or replacement of the Levenhuk product in any country where a Levenhuk office is located if all the warranty conditions are met.

For further details, please visit: [levenhuk.com/warranty](https://www.levenhuk.com/warranty)

If warranty problems arise, or if you need assistance in using your product, contact the local Levenhuk branch.

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Levenhuk® is registered trademark of Levenhuk, Inc.
© 2006–2025 Levenhuk, Inc. All rights reserved.
Levenhuk Inc. (USA): 928 E 124th Ave. Ste D, Tampa, FL 33612, USA, +1 813 468-3001, contact_us@levenhuk.com