



# EBC-X Battery Tester (User Manual V1.6.5)

Thank you for selecting the EBC-X Battery Tester. The battery testers of this series are designed for battery testing and cell grading and can be adopted to a wide range of batteries like Lead-Acid, Li-Fe and Li-Ion.

All devices are shipped fully tested and calibrated and, with proper use, will provide years of reliable service. Please read this user's manual for proper use of the testing system. For any update please visit www.zketech.com to learn more information.

Release Date 2016-12-15 Copyright (C) 2006-2016 ZKETECH



## Contents

1. Features	3
2. Specifications	4
3. Test connection	5
4. Software installation	5
5. Software operation	6
6. Other functions	3
7. Precautions	5



## 1. Features

#### 1.1 Range of application

The EBC-X battery testers can be used for testing and grading a wide range of batteries like LiFe, LiPo, Pb and etc. We provide testers of three categories: high power testers, which cover a range of 5V20A-200A and allow you to test batteries with high capacities, power batteries and super capacitors; standard testers with a range of 5V5A-10A, and low power testers of 5V1mA-1000mA including button cells and other batteries with low capacities. We also provide customaries testers to meet your special needs.

#### 1.2 Hardware characteristics

The testing system consists of upper computer (PC) and lower computer (control circuit unit) and ensure reliable tests.

Each channel works independently and parameters can be set individually in charging and discharging.

The system adopts modularized hardwires which are easy to extend and maintain.

The EBC-X testers connect with independent clamps which can be replaced in testing different batteries. Power and channel indicators on testers allow you to read working state easily.

#### 1.3 Software characteristics

In the software you can read working state, voltage, current, capacity, energy, cycle number and etc. It allows you to set parameters up to 10 steps and for 1000 cycles. Each channel separates from the others and you can select single channel or multi channels. Charging and discharging current, voltage and capacity display in curves which supports for later analization. The software supports automatic data recording and parameter calibration.



# 2. Specifications

Model	X0501	X0510			
Voltage Range	0-4.500V				
Voltage Accuracy		$0.1\%\pm0.003\textrm{V}$			
Current Range	0.001-1.000A		0.1-10.00A		
Current Accuracy	0.2%±0.0005A		0.2%±0.005A		
Charge Current	0.001-1.000A	0.001-1.000A 0.1-5			
Charge Mode	CC-CV				
Charge Termination	Time/Cut-off Current/Manual				
Discharge Current	0.001-1.000A		0.1-10.00A		
Discharge Mode		CC	6		
<b>Discharge Termination</b>	Time	/Cut-off Voltage/Ma	anual		
Testing Channel		8 separate channels			
Heat Dissipation	Fan				
Power Supply	200-240V AC				
Packed Size	300*260*110	22.00	350*260*110		
Packed Weight	3Kg	1 SILV	5.5Kg		



## 3. Test connection

1) Front panel



Each channel has a yellow LED indicator and a 4-core aviation plug socket. The LED indicator lights up during charging or discharging and goes off under other conditions. V+ and V- (as showed in the diagram above) allow you to test voltage while A+ and A- for current testing

2) Back panel

On the back panel you will find power switch (3-core AC) on the right and USB port (4-core aviation plug socket) on the left.

## 4. Software installation

- 1) System requirements
- <sup>2</sup> Processor main frequency>1G Hertz
- **2** RAM>1G, more RAM is suggested.
- **2** Free disc space>1G
- <sup>2</sup> Either Windows XP/Win7/Win8/Win10
- <sup>2</sup> Computer monitor: 1024×768 or higher
- 2) Software installation

First install the tester software. Please double click this icon to start installation.



50)	

Next install the USB port driver. Please double click this icon to start installation.



After installation completes, plug in the USB cable and you will find corresponding port number in Device Manager in your computer (EB software can't identify any port number larger than COM 10, please change the port number to smaller than 10 if it's larger than 10).

## 5. Software operation

1) Switch on the tester.

2) Plug in the USB cable.

3) Connect the tester with a computer: run EBX software and click on Connect Device in the drop-down list of Tools.

File	Tools Help		
	Connect Device		5
	Display Setting		
	Data Analysis		
	Auto Save		
		0.000 An 0.000 An 0.000 An 0.000 An 0.000 An 0.000 An	
		0.000 V	
		0.000 A 0.000 A 0.000 A 0.000 A 0.000 A 0.000 A	
		0.000 Ah	
		[[0.000 Wh][0.000 Wh][0.000 Wh][0.000 Wh][0.000 Wh][0.000 Wh]	

4) In the interface of "Connect to Devices", select the right port (same COM number as that in the computer Device Management) and click on Connect to build connection between the tester and the computer. It allows you to select device numbers in "Address" to include up to six eight-channel testers working through this single COM at the same time.





5) As soon as connection builds, you will see the main interface. In the upper left of the interface you will find a Device List with information of every tester in use.

2016 Tools Help			n
Device List PCV DeviceOff FDC (510 (47V) DeviceOff EBC -(510 (43V)	11 Stp         12 Stp         12 Stp         12 Stp         12 Stp           0.000 A         0.000 Y         0.000 Y         0.000 Y         0.000 A           0.000 A         0.000 A         0.000 A         0.000 A         0.000 A           0.000 Ah         0.000 Ah         0.000 Ah         0.000 Ah         0.000 Ah	4 Sh +         1 + Sh +         1 + Sh +         1 + Sh +         0 + Sh +         <	Y A A A h
	0.000 wh         0.000 wh         0.000 wh         0.000 wh         0.000 wh           24 5/3p         25 5/3p         25 5/3p         25 5/3p         25 5/3p         26 5/3p           0.000 wh         0.000 wh         0.000 wh         0.000 wh         0.000 wh         0.000 wh           0.000 wh         0.000 wh         0.000 wh         0.000 wh         0.000 wh         0.000 wh           0.000 wh         0.000 wh         0.000 wh         0.000 wh         0.000 wh         0.000 wh	1 Geo         25 Step         25 Step         0.000 Wh	Coc V A A Wh

6) Fix batteries on the testing clamps, and voltage of each channel will display in the right part of the interface.



7) Right click on any channel to set parameters in Parameter Setting.

1-1 Stop 3.049 V 9 9 1-2 Stop 3.049 V 9 9 9 9 9 9 9 9 1-2 Stop 3.049 V 9 9 9 9 9 9 9 9 9 9 9 9 9	1-3 Stop 3.049 V 0.000 A 0.000 Ah 0.000 Wh	1-4 Stop 3.049 V 0.000 A 0.000 Ah 0.000 Wh	<ul> <li>1-5 Stop</li> <li>4.600 V.</li> <li>0.000 A</li> <li>Volt Ah</li> <li>Error Wh</li> </ul>	1-6 Stop 3.049 V 0.000 A 0.000 Ah 0.000 Wh	1-7 Stop 3.049 V 0.000 A 0.000 Ah 0.000 Wh	1-8 Stop 3.049 V 0.000 A 0.000 Ah 0.000 Wh
Channel Continue Channel Save Device Start Device Stop Device Continue	2-3 Stop 3.107 V 0.000 A 0.000 Ah 0.000 Wh	2-4 Stop 3.107 V 0.000 A 0.000 Ah 0.000 Wh	<ul> <li>2-5 Stop</li> <li>4.600 V</li> <li>0.000 A</li> <li>Volt Ah</li> <li>Error Wh</li> </ul>	2-6 Stop 3.107 V 0.000 A 0.000 Ah 0.000 Wh	2-7 Stop 3.107 V 0.000 A 0.000 Ah 0.000 Wh	2-8 Stop 3.107 V 0.000 A 0.000 Ah 0.000 Wh

8) Set parameters in the Parameter Setting interface. It allows you to set Step-Mode-Current-Voltage and etc.



<	Parameter Setting						
	Step	Mode	Current	Voltage	Cut-off	Time/Cir	
	1	C-CV	02.00	04.20	00.10	0	
	2	Wait	00.00	00.00	00.00	5	
	3	D-CC	02.00	02.80	00.00	0	
	4	Wait	00.00	00.00	00.00	5	
	5	Cycle	00.00	00.00	00.00	3	
	6						
	7						
	8						
	9						
	10						
	Step: Step6 - Mode: Cycle - Cycle 1						
	Curr: 00.00 A Volt: 00.00 V Cut-off: 00.00 A						
	Clear Set Open Save Copy						
	Start S	itep: 1	Channel:	1	• ОК	Cancel	

- <sup>2</sup> Clear: Remove all parameters you've set.
- <sup>2</sup> Set: Set test parameters for the current step.
- <sup>2</sup> Open: Open saved files of parameters.
- <sup>2</sup> Save: Save current parameter files to local disks.
- <sup>2</sup> Copy: Copy parameters of other steps to current step.
- <sup>2</sup> Start Step: Select the starting step.
- <sup>2</sup> Channel: Select the channel(s) to use the current parameters.



**2** OK: Use parameters to selected channel.



<sup>2</sup> Cancel: Cancel parameter setting.

9) Right click on a channel and select "Channel Start" on the pop-up menu to test that channel, or select "Device Start" to test all the eight channels in a tester.



10) Double click on a channel to open the interface with detailed data of that channel.



Double click on the curve display area to go back to the main interface to have an overall view of parameters in every channel. The colorful box on the top of each channel means as follows:



- 1-1 C-CV	= 1-2 C-CV	= 1-3 C-CV	= 1-4 C-CV	= 1.5 C-CV	= 1-6 C-CV	= 1.7 C-CV	= 1-8 C-CV
3.029 V							
0.541 Ah	0.767 Ah	0.523 Ah	0.530 Ah	0.646 Ah	0.632 Ah	0.593 Ah	0.441 Ah
3.338 Wh	3.102 Wh	3.351 Wh	3.214 Wh	3.283 Wh	2.905 Wh	3.378 Wh	3.301 Wh

**2** Green **=**: Discharging

- **2** Red **=**: Charging
- **2** Yellow **=**: Waiting
- 2 Gray 🔳: Stop

The system allows you to change fonts and colors of parameters displayed on the right of the main interface. It takes a single click on characters at the bottom.

Volt:	Curr:	Can:	Ene:
VOIL	Curr.	🖡 Cap. 💶 🔤	

11) You can switch between "Graph View" and "List View" by clicking on the icons.



12) After testing completes, you can identify channels with required capacities by clicking on "Cap Filter" to open the menu of "Capacity Search".



Set capacity ranges and click on "OK" to display certain channels and click on "Clear" to return to original displays.



13) Click on "Data Report" at the bottom of the main interface and you will see all test results.

<u></u>	- Alaport									
Channel	Step	Capacity	Step	Capacity	Step	Capacity	Step	Capacity	Step	Capa_
1-1	1-1 C-CV	1.582								
1-2	1-1 C-CV	2.039								
1-3	1-1 C-CV	1.792								
1-4	1-1 C-CV	1.82								
1-5	1-1 C-CV	2.01								
1-6	1-1 C-CV	1.879								
1-7	1-1 C-CV	1.888								
1-8	1-1 C-CV	1.864								
										2.14
						Ta:		10		
21=1										
<u>.                                    </u>										

14) In the diplay of "Graph View", you can right click on a certain channel to open the drop-down menu. On the menu, select "Channel Save" or "Device Save" to save data to hard disks.



- 15) Options on the drop-down menu
- <sup>2</sup> Channel Start: start testing the present channel
- <sup>2</sup> Channel Stop: stop testing the present channel
- <sup>2</sup> Channel Continue: continue last test in the present channel
- <sup>2</sup> Channel Save: save testing results of the present channel
- <sup>2</sup> Device Start: start testing all channels
- <sup>2</sup> Device Stop: stop testing all channels



- <sup>2</sup> Device Continue: continue testing in all channels of the device
- <sup>2</sup> Device Save: save testing results of every channel in the device
- <sup>2</sup> Multiple Start: start testing selected channels
- <sup>2</sup> Multiple Stop: stop testing selected channels
- <sup>2</sup> Multiple Continue: continue last testing in selected channels
- <sup>2</sup> Parameter Calibration: calibrate voltage and current of the present channel
- 16) Parameter Calibration

🏶 CalibrationDevicel Channell 🔀						
Voltage Calibration						
Low Volt:	0.000	۷	Calibration			
High Volt:	0.000	۷	Calibration			
Current Calibration						
Low Curr:	0.000	А	Calibration			
High Curr:	0.000	А	Calibration			
OK						

To calibrate the voltage, you need a high precision multimeter and an adjustable voltage source (can be replaced by two batteries of different voltage).

Connect the channel to a low voltage. Measure the actual value with the multimeter (for example 1V) and fill it in "Low Volt". Click "Calibration".

Then connect the channel to a high voltage. Measure the actual value with the multimeter (for example 4V) and fill it in "High Volt". Click "Calibration".

Click "OK" to save the voltage calibration parameters. A small difference between the high and the low voltage is enough.

To calibrate the current, you need a high precision multimeter and a battery.

Connect the battery to the channel, set a low current and start discharging it. Measure the actual value with the multimeter (for example 0.5A) and fill it in "Low Curr". Click "Calibration".

Set a high current and start discharging. Measure the actual value (for example 2.5A) and fill it in "High Curr". Click "Calibration".

Click "OK" button to save discharging currents. A small current difference is enough.



Follow the same steps and complete calibrating charging current. Click "Cancel" button to go back to the main interface.

# 6. Other functions

1) Auto Save



Click on Tools to open the drop-down menu and you can select Auto Save to save testing results automatically in the Dada file in the installation directory.

2) Display Setting

🥥 Display Sett	ing 🔀
Color	Line
Grid:	Majoy: Solid 👻
Bq:	Minor: Dot 💌
Axes:	Line: 1 🔹
Curr:	-l egend
Volt:	
(	OK Cancel
	OK Cancel

Click on Tools to open the drop-down menu and you can select Display Setting to set Curves and Lines in curve display area.

## 3) Data Analysis

Click on Tools to open the drop-down menu and you can select Data Analysis to open the menu of Curve Data. On the menu you can select data files by clicking on File. The default display of files contains testing data and curves of each testing step. You can select curve types at the right bottom and click on Refresh to get curves you want.

🗘 Cu	urve Data	
File	Curve	
Ope	en Data	
Sat	ve Bata	



#### **EBC-X Battery Tester**

Cwr (V) 50 gr	**																			Second Second
M) 50 m																				
M) 50 m					- 2	EDO.	VARIA	-						3255	Oycle	Step	Mode	Capacity	Energy	Time
50 J.						EBU	×0510	1-1						A	1	1	C-CV	2.366	9.323	80.22.1
10.0	100 100	. itin	dine.		1.1	1.36	10		0000	1.00	and an		1.10	10.50	1	2	Wait	0.000	0 000	00.05.0
m .[[	The second secon					100							CITE I C	9.45	1	3	D-CC	2.363	8.949	00:14:1
~][		- 51	den de			1.02				1 3				1.000	1	4	Wait	0.000	0.000	00.05.0
10 -	- 1				1	1			-1	1 3			-1	8.40	2	1	C-CV	2.363	9.311	00:22:00
. II	-/		\	14	1	1.64	<b>\</b>		-1		A		1	10000	2	2	Wait	0.000	0.000	00.05:0
0 -	1	- 1	1			1	1		1	3	1		1	- 7.35	2	3	D-CC	2.363	8.950	00.14.1
			$\mathbf{X}$	1	111	( Comp		11/	1.1		$+ \lambda -$		and the	1000	2	4	Wait	0.000	0.000	00.05.0
0 1										1				- 6.30	3	1	C-CV	2.361	9.303	00.22.1
	1.1	100		11	1251	173	1720		115	1	111 22		1.51	100	3	2	Wait	0.000	0.000	00:05:0
~ ][	and and	1.2	1.1		1.97	1.11	a la contra		11.2	1.1.1				249	3	3	D-CC	2,360	8,939	00.14.1
					1									- 120	3	4	Wait	0.000	0.000	00.05.0
~ <b> </b>		and a	-			and .			orter				- aler	1222	4	1	C-CV	2.339	9.218	00 17 4
0 -					-	1			- 11					- 3.15						
	and and			(	-			1	1.222			¥1		1		_	-			
90 -					1			1				1	- 5	- 210						
	1		1000	1	100	1000			1997	1	1	t t		1						
70.1	100				160				0.00		1.000			1.05						
, li	13	V	1.1			X			1	N				Lam						
00.00	200 0015 M	5.41 00 ode	31.23 (	0.47.0 Cu	75 OT	02.46 (	01:18:28 Volt	01.34 age	1.09 01	.49:50 O	02:05:32 ut-off	02.2	1.13 02 Tim	36.55 8/Oyc						
CTT	D	-CC		10	0.00		04	20		0	0.10			0						
2	Ŵ	Vait		00	00.0		00	00		0	0.00			5						
3	0	CV		10	0.00		02	70		0	0.00			0						
1	W	Vait		00	0.00		00	00		0	0.00			5						
5	9	ycle		.00	0.00		00	00		0	0.00			5	-					
														0	C DS	C-Tim	e C CH	G-Time @	ALL-Time	Refre



WWW.ZKETECH.COM Copyright (C) 2006-2016 ZKETECH 1

in the



	Data										
e Curv	**										
M		EBC-X05	10 1-1 CHG-Tin	ne	ivi	Cycle 1	Step 1	Mode C-CV	Capacity 2.366	Energy 9.323	Time 00.22.14
140 7	edine has been	And the loss	in hack of a	and the basis	- 10.50	2	1	C-CV	2.363	9.311	00.22.02
125					9.45	3	1	C-CV	2.361	9.303	00.22:19
			anness anness			4	1	O-CV	2.339	9.218	00.17:47
10 -	1 1 1	+ + +			- 8.40						
~ <b>H</b>	and and and a second	apresident in the			an and a second second						
95					- 7.35	-					
		50	100 C		6.20	-					
~1	1	dan dan da			- 0.30	-	-	<u>1</u>			
65 -	1				- 5.25						
	1.1	den minimi		an a							
50 -		2 2 2	1	1 1 1	- 4.20						
. II	. Jacobardian des	1000	1	이 이번 이 이번							
21	State and the second				- 3.15						
20 -					- 210						
	Ob. 11. (10)			A 10 10 10 10 10 10 10 10 10 10 10 10 10							
6526	and the second			Northern confe							
.05 -				X	+ 1.05						
05 - 90 -				X	- 1.05						
05 90 00.00	100 00 02 14 00 04 2 Mode	7 00:06:41 00:08:55 Current	001106 001322 ( Voltage	00.15.36 00.17.50 0 Cut-off	1.05 0.20.03 00.22.17 Time/Cyc						
05 90 00.00	100 0010214 001042 Mode D-CC	7 00.06.41 00.08.55 Current 10.00	0011108 001322 0 Voltage 0420	00.15.36 00.17.50 0 Cut-off 00.10	- 1.05 0.20.03 00.22.17 Time/Cyc 0						
00.00 30 50000 51ep 1 2	100 00 02 14 00 04 2 Mode D-CC Wet	7 00.06.41 00.08.95 Current 10.00 00.00	0011.08 00.13.22 ( Voltage 0.4.20 00.00	00.15.36 00.17.50 0 Cut-off 00.10 00.00	1.05 0.00 00 22:17 Time/Cyc 0 5						
105 290 00.00 Step 1 2 3	100 00 02 14 00 04 2 Mode D-CC West C-CV	7 00.06.41 00.08.95 Current 10.00 00.00 18.00	0011.08 00.13.22 0 Voltage 04.20 00.00 02.70	00.15.36 00.17.50 0 Cut-off 00.10 00.00 00.00	1.05 0.000 0.2003 00.2217 Time/Cyc 0 5 0						
00.00 290 305 290 290 5tep 1 2 3 4	Mode D-CC West C-CV Wat	7 00.06.41 00.0855 Current 10.00 00.00 10.00 00.00	0011.08 00.13.22 0 Voltage 04.20 00.00 02.70 00.00	00.15.36 00.17.50 0 Cut-off 00.10 00.00 00.00 00.00	1.05 0.000 00.2217 Time/Cyc 0 5 0 5 0 5						
3.05 2.90 00.00 Step 1 2 3 4 5	Mode D-OC Weat C-CV Wat Cycle	7 00.06.41 00.08.55 Current 10.00 00.00 18.00 00.00 00.00 00.00	0011.08 0013.22 0 Voltage 04.20 00.00 02.70 00.00 00.00	00.15.36 00.17.50 0 Cut-off 00.10 00.00 00.00 00.00 00.00	1.05 0.20.03 00.22.17 Time/Cyc 0 5 0 5 5	0					
00.00 300p 1 2 3 4 5	Mode D-CC Weat C-CV Wat Cycle	7 00.06.41 00.08.95 Current 10.00 00.00 10.00 00.00 00.00 00.00	0011.08 00:13:22 0 Voltage 04:20 00:00 02:70 00:00 00:00	00.15.36 00.17.50 0 Cut-off 00.10 00.00 00.00 00.00 00.00 00.00	1.05 0.20103 00:22:17 Time/Oyc 0 5 5 5 5	C DS	cC-Tim	е « сн	G-Time f	ALL-Time	Bafach

## 7. Precautions

1) Please switch on the tester before you connect it to your computer.

- 2) Please remove all batteries from your tester before you switch it off.
- 3) Not support testing batteries in series.
- 4) Always use the tester within the allowed range.
- 5) Follow calibration instructions in case unproper operation reduce precision.

We always focus on innovation and improvement and will keep upgrading the software. Please frequently visit our official website: www.zketech.com to download the latest manuals.