F-1501 Wanda Plaza, No.18 Shijingshan Road, Beijing 100043, China Tel:+8610 88696020 Fax:+8610 88696085

A33.5006

LCD Digital Microscope

Catalogue



Please read this manual carefully before use microscope!



PRODUCT DESCRIPTION

6.6				
A33.5006 4.3-inch LCD Inspection Microscope Specification				
Magnification Ratio	600x			
Still Capture Resolution	8M, 5M, 3M, 1.3M, VGA			
Focus Range	Manual Focus From 5mm To 80mm			
Frame Rate	Max 30f/s Under 600 Lus Brightness			
Video Format	AVI			
Photo Format	JPEG			
Light Source	LEDs (Brightness Adjustable)			
Battery	Built-in Li-ion Battery			
DC Input	5V/1A			
TV-out	To Any Monitor With TV-in			
Storage Micro-SD Card				
Bundle Software	Software Portable Capture Plus With Measurement			
Operation System For Bundle Software	Windows 7/8/10			
	Mac OS X 10.12 Or Above			
Software	Android Smartphone/Tablet			

www.cnoec.com.cn Page 2 of 7 sale@cnoec.com


Product Size	116mm (L) x 95mm (W)	
Weight	Microscope-220g; Stand-220g	
	Microscope	
	Observation Stand	
Doolsono	Software CD	
Package	USB Cable	
Contents	TV Cable	
	Micro-SD Card (1GB)	
	User's Manual	







PART NAME AND FUNCTION




1	Focus Knob	9	Micro-SD Card Slot
2	Menu	10	Digital Transform Multiples
3	Video, Photo Capture, Playback Switch	11	Battery Charging Indicator Light (Red For Charging,
4	Up		Full Of Rear Light Out)
5	Down	12	Television Outlet
6	Video & Photo Capture	13	Charging Port: Charge The Internal Battery
7	Boot to Turn It Off		USB Connected Computer: Read The Card, Or
8	LED Source Adjust		Enable Live Video On Your PC/Android Device

APPLICATION AND CASE:

The 4.3-inch LCD Standalone Inspection Digital Microscope is a fantastic microscope that provides a good way for you to share high definition microscopic observation with others. The 4.3" quality LCD panel brings vivid crystal image.

Equipped with professional microscopic lens, it helps you to see tiny objects easily. The magnification can be up to 600x. Built-in rechargeable Li-ion battery provides a 6-7rs standalone working time. The 8 bright LEDs, with brightness adjustable, are bright enough for perfect observation.





PACKING AND SHIPPING

Packing: Microscope-220g; Stand-220g;

Product Size 116mm (L) x 95mm (W)

Shipping Time: 3-10 Days From Stock

A33.5006



















We Remain Right To Change Technical Data Without Notice Due To Upgrade.